




Model No.




Example:

QHF-SV06-X-T6-K1.2-S refers to SV06 6 port selector valve with orifice 1.2mm and sapphire valve spool


## Technical Parameter

|  |   |               |               |
|--|---|---------------|---------------|
| Valve configuration                              | 6/8/10 port   | 12 port       | 16 port       |
| Orifice  | 1.2mm   | 1.0mm         | 1.0mm         |
| Fluid path                                       | PCTFE valve head + Sapphire valve spool   |               |               |
| Dead volume                                      | 5.41 $\mu$ L  | 6.08 $\mu$ L  | 10.4 $\mu$ L  |
| Port to port volume                              | 27.5 $\mu$ L  | 22.43 $\mu$ L | 33.68 $\mu$ L |
| Max drive power (torque)                         | 2.6N/m  | 3.0N/m        | 3.5N/m        |
| Secondary drive power (torque)                   | 0.32N/m   | 0.36N/m       | 0.42N/m       |
| Pressure rating                                  | 0-1.0Mpa (air) 0-1.6Mpa (water)   |               |               |
| Detection  | auto-detect original position when powered on<br>(when selectively open or close the valve) |               |               |
| Liquid temperature                               | 0-150°C   |               |               |
| Connection                                       | 1/4-28UNF Female  |               |               |
| Replaceable parts                                | stator replaceable, sealed rotor unreplaceable  |               |               |
| Transposition                                    | random start to different ports   |               |               |
| Driver   | non-optional (Integrated)   |               |               |
| Switch speed                                     | $\leq$ 5s/circle  |               |               |
| Communication                                    | RS232/RS485: 9600bps, 19200dps, 38400dps, 57600dps, 115200dps                               |               |               |
|  | CAN: 100Kbps, 200Kbps, 500Kbps, 1Mbps   |               |               |
| Address & Parameter settings                     | Communication Interface   |               |               |
| Power supply                                     | DC24V/3A  |               |               |
| Maximum power                                    | 60W   |               |               |
| Working environment                              | -10°C - 50°C  |               |               |
|  | $\leq$ 80% relative humidity, non-condensing  |               |               |
| Dimension (L*W*H)                                | 60*51*145.8mm   | 60*51*156.5mm | 60*51*175.5mm |
| Net weight (kg)                                  | 0.73kg  | 0.86kg        | 1.02kg        |
| OEM support for special application requirements |   |               |               |


### Port to Port Volume



6/8/10 Port





12 Port



16 Port

### Dead Volume



## Valve Configuration



6 Port

8 Port

10 Port


12 Port

16 Port

## Valve Reset

**Reset Direction:** Counterclockwise

**Reset Position:** Between port 1 and maximum port number (6/8/10/12/16), when rotor at reset position, it does NOT connect to any other ports.







## Dimension (unit: mm)




SV-06 Dimension.pdf

## Mounting Size (unit: mm)



Please contact customer service if you need 2D and 3D file, thanks for cooperation!

## Driver Port



| Port Name | Function               |
|-----------|------------------------|
| +         | DC24V positive         |
| -         | DC24V negative         |
| TX        | RS232 TX               |
| RX        | RS232 RX (Data output) |
| GND       | GND                    |
| H         | CANH                   |
| L         | CANL                   |
| A         | RS485 A                |
| B         | RS485 B                |

## Driver Control Instruction

The data between selector valve and upper monitor (PC, PLC, Raspberry Pi, micro-controller) was transmitted by serial communication RS232/RS485/CAN

Communication Form: Asynchronous serial communication; Command and data frames are sum check 2 Byte; Commands and data are hexadecimal numbers; Command parameters saved by little-endian mode.

Communication Interface: RS232 or RS485 or CAN

Communication Mode: Bidirectional asynchronous; master-slave mode

Baud rate: 9600bps,19200bps,38400bps,57600bps,115200bps (RS232/RS485) / 100K, 200K, 500K, 1M (CAN)

Data bit: 8

Even-odd Check: None

Response Time: <1 second

**Common Command (send command 8 bytes / received command 8 bytes)**

### Send Command (Monitor)

| B0         | B1      | B2       | B3               | B4       | B5       | B6        | B7        |
|------------|---------|----------|------------------|----------|----------|-----------|-----------|
| Start code | Address | Function | Status Parameter |          | End code | Sum check |           |
| STX        | ADDR    | FUNC     | 1-8 bit          | 9-16 bit | ETX      | Low byte  | High byte |

First byte "STX" = Start code (CCH)

Second byte "ADDR" = Address of slave device (0x00~0xFF)

Third byte "FUNC" = Function code

Forth Fifth byte = Status parameters of function code

Sixth byte = End code (DDH)

Seventh Eighth byte = Sum check code from byte 1 to 6

## Response Command (Slave)

| B0         | B1      | B2       | B3               | B4       | B5       | B6        | B7        |
|------------|---------|----------|------------------|----------|----------|-----------|-----------|
| Start code | Address | Function | Status Parameter |          | End code | Sum check |           |
| STX        | ADDR    | FUNC     | 1-8 bit          | 9-16 bit | ETX      | Low byte  | High byte |

First byte "STX" = Start code (CCH)

Second byte "ADDR" = Address of slave device (0x00~0xFF)

Third byte "FUNC" = Function code

Forth Fifth byte = Status parameters of function code

Sixth byte = End code (DDH)


Seventh Eighth byte = Sum check code from byte 1 to 6

## Common Command

| Command Type    | Command code | Instruction  | Parameter (B3, B4)     | Response Parameter (B3, B4)   |
|-----------------|--------------|--|------------------------|---|
| Query Command   | 0x21         | Query RS232 baud rate                              | 0x0000                 | Factory default 9600bps (B4 =0x00)<br>B3=0x00 baud rate 9600bps<br>B3=0x01 baud rate 19200bps |
|                 | 0x22         | Query RS485 baud rate                              | 0x0000                 | B3=0x02 baud rate 38400bps<br>B3=0x03 baud rate 57600bps<br>B3=0x04 baud rate 115200bps       |
|                 | 0x23         | Query CAN baud rate                                | 0x0000                 | B3B4=0x0000 100Kbps<br>B3B4=0x0001 200Kbps<br>B3B4=0x0002 500Kbps<br>B3B4=0x0003 1Mbps        |
|                 | 0x2e         | Query automatic reset when power on                | 0x0000                 | (B4=0x00)<br>B3=0x00 Non-auto reset when power on<br>B3=0x01 automatic reset when power on    |
|                 | 0x30         | Query CAN destination address                      | 0x0000                 | B3=0xXX (B4=0x00)<br>XX = 00~FF, default as 00  |
|                 | 0x3e         | Query current located port                         | 0x0000                 | B3B4 = Current position of coded disc   |
|                 | 0x3f         | Query current firmware version                     | 0x0000                 | B3B4 = Software version number (hexadecimal)  |
| Control Command | 0x44         | Motor runs by coded disc, auto-select optimal path | Port No. (1-Max. port) | RS232: B2=0x00<br>B3B4= Internal data (randomly appears)<br>RS485: B2=0xFE B3B4=0x0000        |
|                 | 0x45         | Reset  | 0x0000                 | RS232: B2=0x00<br>B3B4= Internal data (randomly appears)<br>RS485: B2=0xFE B3B4=0x0000        |
|                 | 0x49         | Strong stop  | 0x0000                 | B3B4=Rest steps   |
|                 | 0x4a         | Query motor status                                 | 0x0000                 | B3B4=0x0000   |

| Response Status | Response Code (B2) | Parameter Instruction |
|-----------------|--------------------|-----------------------|
|                 | 0x00               | Normal status         |
|                 | 0x01               | Frame error           |
|                 | 0x02               | Parameter error       |
|                 | 0x03               | Optocoupler error     |
|                 | 0x04               | Motor busy            |
|                 | 0x05               | Motor stalling        |
|                 | 0x06               | Unknown position      |
|                 | 0xfe               | Task suspension       |
|                 | 0xff               | Unknown error         |

**Factory Command (send command 14 bytes / received command 8 bytes)**



**Send Command (Monitor)**

| B0                | B1             | B2              | B3,B4,B5,B6     | B7                        | B8       | B9        | B10             | B11              | B12      | B13       |
|-------------------|----------------|-----------------|-----------------|---------------------------|----------|-----------|-----------------|------------------|----------|-----------|
| <b>Start code</b> | <b>Address</b> | <b>Function</b> | <b>Password</b> | <b>Function Parameter</b> |          |           | <b>End code</b> | <b>Sum check</b> |          |           |
| STX               | ADDR           | FUNC            |                 | 1-8 bit                   | 9-16 bit | 17-24 bit | 25-32 bit       | ETX              | Low byte | High byte |

First byte "STX" = Start code (CCH)

Second byte "ADDR" = Address of slave device (0x00~0xFF)

Third byte "FUNC" = Function code

Forth to Seventh byte = Password of factory command

Eighth to Eleventh byte = Parameters of function code

Twelfth byte = End code (DDH)

Thirteenth to fourteenth byte = Sum check code from byte 1 to 12

### Response Command (Slave)

| B0         | B1      | B2       | B3               | B4       | B5       | B6        | B7        |
|------------|---------|----------|------------------|----------|----------|-----------|-----------|
| Start code | Address | Function | Status Parameter |          | End code | Sum check |           |
| STX        | ADDR    | FUNC     | 1-8 bit          | 9-16 bit | ETX      | Low byte  | High byte |

First byte "STX" = Start code (CCH)

Second byte "ADDR" = Address of slave device (0x00~0xFF)

Third byte "FUNC" = Function code

Forth Fifth byte = Status parameters of function code

Sixth byte = End code (DDH)

Seventh Eighth byte = Sum check code from byte 1 to 6

### Factory Command


| Command Type    | Command code | Instruction                       | Parameter (B3, B4) | Response Parameter (B3, B4)  |
|-----------------|--------------|-----------------------------------|--------------------|--|
| Factory Command | 0x00         | Set device address                | 0x0000~0x00FF      | Address: 0~255   |
|                 | 0x01         | Set RS232 baud rate               | 0x0000~0x0004      | Factory default 9600bps (B4=0x00)  |
|                 | 0x02         | Set RS485 baud rate               | 0x0000~0x0004      | B3=0x00 baud rate 9600bps<br>B3=0x01 baud rate 19200bps<br>B3=0x02 baud rate 38400bps<br>B3=0x03 baud rate 57600bps<br>B3=0x04 baud rate 115200bps |
|                 | 0x03         | Set CAN baud rate                 | 0x0000~0x0003      | B3B4=0x0000 100Kbps<br>B3B4=0x0001 200Kbps<br>B3B4=0x0002 500Kbps<br>B3B4=0x0003 1Mbps   |
|                 | 0x0E         | Set automatic reset when power on | 0x0000~0x0001      | B3B4=0x0000<br>Non-auto reset when power on<br><br>B3B4=0x0001<br>Automatic reset when power on  |
|                 | 0x10         | Set CAN destination address       | 0x0000~0x00FF      | Address: 0~255   |

#### Instructions:

(1) Code B2 in response command means current motor status. Only when B2=0x00 motor works normally.

Other status parameters are in above table.

When device controlled by RS485 and send control command B2-0x44 or 0x45, status parameter in response command is FE (task suspending) which means motor is working as command told, if now send other commands (except for query command), the status parameter in response command will be 04 (motor busy), if send polling command 0x4a, the status parameter in response command will be 00 (motor under normal status).




- ① Send control command (B2=0x44 or 0x45), received code FE (task suspending)
- ② Send query command 0X4A, received code 04 (motor busy)
- ③ Send query command 0X4A, received code 00 (motor under normal status)

(2) Parameter B3,B4 in response command make sense only when send query command; when send setting or control command, response parameters make no sense with default 00 00. When send query command and parameter B2 in response command is 00, then response parameter B3, B4 make sense, received value is the query result. E.g. when send inquiry command 0x21 (query RS232 baud rate), response command B3 B4 = 04 00, it means baud rate of RS232 is 115200bps.

**Note:** all command parameters set and saved by little-endian mode. Little-endian means lower data bit saved in the lower address bit, higher data bit saved in the higher address bit.

### Application Notice

- ◆ Please ensure input voltage was as required
- ◆ Please use original serial port lines
- ◆ Communication RS232, RS485, CAN are under Non-isolation mode, hot swapping unsupported.
- ◆ Please cover the unused ports with suitable coned plugs when laid aside to avoid impurity substance and air
- ◆ Please don't depart all parts on the device and keep all the labels safe and sound in case of warranty service
- ◆ Please read above operation instructions and communication protocols carefully, do not input data randomly.
- ◆ When control several instruments through RS485 and CAN, impedance matching need to be considered.



Nanjing Runze Fluid Control Equipment Co.,Ltd  
 No.9 Tianxing West Road Dongshan Street Jiangning District  
 Nanjing City, Jiangsu Province, China  
 Mobile: +86 17366384502  
 Email: [min.zhu@runzeliuti.com](mailto:min.zhu@runzeliuti.com)  
 Contact: Julie Zhu