


# Color Flangeless Fittings


## Features:

- Different fitting colors for easy differentiation of flow paths
- Imported PP material, non-toxic, odorless, heat resistant, corrosion resistant, widely used in biological medicine, environmental protection, food industry
- Standard threads and curvature, no gap, no wobble and good sealing with the corresponding threads
- Available in five sizes of 1.6/2.0/2.5/3.0/3.2mm OD rigid tube, with a pressure capacity of 7 kg
- Quick and easy installation, no tools required, cost-effective



### ★ Application Note


Installation steps same as flangeless fittings




※ Make sure the ferrule is oriented correctly! The narrow end of the ferrule should towards the nut.

### Nut

Unit: mm



### Ferrule



**YC-016**  
For 1.6mm OD tube  
**YC-020**  
For 2.0mm OD tube

**YC-025**  
For 2.5mm OD tube  
**YC-030**  
For 3.0mm OD tube  
**YC-032**  
For 3.2mm OD tube

Model No.	Thread	Tube OD	Material
<b>Nut</b>			
WD10A-1/4-28-P-Red	1/4-28UNF	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PP Red
WD10A-1/4-28-P-Dark Blue	1/4-28UNF	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PP Dark Blue
WD10A-1/4-28-P-Beige Blue	1/4-28UNF	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PP Beige Blue
WD10A-1/4-28-P-Deep Purple	1/4-28UNF	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PP Deep Purple
WD10A-1/4-28-P-Lemon Yellow	1/4-28UNF	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PP Lemon Yellow
WD10A-1/4-28-P-Dark Brown	1/4-28UNF	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PP Dark Brown
WD10A-1/4-28-P-Grass Green	1/4-28UNF	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PP Grass Green
WD10A-1/4-28-P-Deep Black	1/4-28UNF	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PP Deep Black
WD10A-M6-P-Red	M6	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PP Red
WD10A-M6-P-Dark Blue	M6	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PP Dark Blue
WD10A-M6-P-Beige Blue	M6	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PP Beige Blue
WD10A-M6-P-Deep Purple	M6	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PP Deep Purple
WD10A-M6-P-Lemon Yellow	M6	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PP Lemon Yellow
WD10A-M6-P-Dark Brown	M6	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PP Dark Brown
WD10A-M6-P-Grass Green	M6	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PP Grass Green
WD10A-M6-P-Deep Black	M6	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PP Deep Black
WD10A-1/4-28-E-Red	1/4-28UNF	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PTFE Red
WD10A-M6-E-Red	M6	1.6mm/2.0mm/2.5mm/3.0mm/3.2mm OD	PTFE Red
<b>Ferrule</b>			
YC-016-P		1.6mm OD	PP natural
YC-020-P		2.0mm OD	PP natural
YC-025-P		2.5mm OD	PP natural
YC-030-P		3.0mm OD	PP natural
YC-032-P		3.2mm OD	PP natural