


Female Thread Adapters

Features:

- Standard threads and curvature, no gap, no wobble and good sealing with the corresponding threads
- Imported raw material, non-toxic, tasteless, heat and corrosion resistant.
- PP material is widely used in biomedical, environmental protection and food industries; PTFE material is also suitable for chemical industry and other fields




Application Note




※ Can be connected to 1/4-28 UNF male threaded fittings to allow for any connection between rigid tube and rigid tube, rigid tube and hose, hose and hose of the same or different tube diameters


Unit: mm




NF08-28P/E
Female Union



NY03-28P/E
Female Y Adapter



NT03-28P/E
Female T Adapter



NS04-28P/E
Female Cross Adapter

Model No.	Thread	Material
Female Union		
NF08-28P	1/4-28UNF female thread	PP natural
NF08-28E	1/4-28UNF female thread	PTFE natural
Female Y Adapter		
NY03-28P	1/4-28UNF female thread	PP natural
NY03-28E	1/4-28UNF female thread	PTFE natural
Female T Adapter		
NT03-28P	1/4-28UNF female thread	PP natural
NT03-28E	1/4-28UNF female thread	PTFE natural
Female Cross Adapter		
NS04-28P	1/4-28UNF female thread	PP natural
NS04-28E	1/4-28UNF female thread	PTFE natural