

Gastight Syringe

气密性进样器



- Standard 1/4-28UNF thread, good tightness
- Borosilicate syringe, PTFE piston, 304 SST pushrod high precision, corrosion resistance
- Easy clean and sterilization, reusable
- Flexible connections to meet different applications
- Maximum pressure rating 0.68Mpa (6.8Bar)


RZ-JYQ - 0.5 - W - 1 - U

Model No. Volume Thread Channel Thread


①	0.5	0.5ml	10	10ml	②	W	Male	③	1	single channel	④	U	1/4-28UNF
	1	1ml	25	25ml		N	Female		2	double channel		M6	M6
	2.5	2.5ml	50	50ml								M10	M10
	5	5ml	100	100ml									

Connection

1. Rigid Tubing



2. Flexible Tubing




3. Needle



Specification (OEM support)

Model No.	Volume	Channel	Thread	Syringe ID	Syringe OD	Stroke	Net weight
RZ-JYQ-0.5-W-1-U	0.5ml	1(single)	Male	5mm	8.1mm	25.5mm	16.4g
RZ-JYQ-1-W-1-U	1ml	1(single)	Male	5mm	8.1mm	50.93mm	21.6g
RZ-JYQ-2.5-W-1-U	2.5ml	1(single)	Male	7.25mm	10.8mm	60.56mm	38g
RZ-JYQ-2.5-N-1-U	2.5ml	1(single)	Female	7.25mm	10.8mm	60.56mm	38g
RZ-JYQ-5-W-1-U	5ml	1(single)	Male	10.3mm	15.18mm	60mm	89g
RZ-JYQ-5-N-1-U	5ml	1(single)	Female	10.3mm	15.18mm	60mm	79g
RZ-JYQ-10-W-1-U	10ml	1(single)	Male	14.55mm	18.6mm	60.14mm	106g
RZ-JYQ-10-N-1-U	10ml	1(single)	Female	14.55mm	18.6mm	60.14mm	85.4g
RZ-JYQ-25-N-1-U	25ml	1(single)	Female	23mm	29mm	60.17mm	183g
RZ-JYQ-25-N-2-U	25ml	2(double)	Female	23mm	29mm	60.17mm	200g
RZ-JYQ-50-N-1-U	50ml	1(single)	Female	32mm	39.4mm	60.17mm	280g
RZ-JYQ-50-N-2-U	50ml	2(double)	Female	32mm	39.4mm	60.17mm	297g
RZ-JYQ-100-N-1-U	100ml	1(single)	Female	34mm	40.2mm	110.14mm	369g
RZ-JYQ-100-N-2-U	100ml	2(double)	Female	34mm	40.2mm	110.14mm	386g



OEM Service

