




Luer Adapters

软管接头

Luer adapters were designed to work in a variety of applications such as medical and laboratory instruments. You can create your own quick connect union or adapter by connecting any male luer to any female luer.

- Quick disconnect and replacement
- Adapts luers to 1/4"-28 threaded ports
- Medical grade PP version (Other material customize support)
- For Flexible tubing ID 1/16" 3/32" 1/8" 5/32" 3/16" 1/4"
- For Rigid tubing OD 1/16" (1.6), 2.0/2.5/3.0/3.2mm, 1/8" (3.175)

Model No.




RH-G032: Male luer adapter suit tubing ID2.4-3.2mm

RH-2M: Female luer-to-luer adapter

RH-M-(1/4-28)W: Female luer-to-male thread 1/4-28UNF

Connection




Specification (unit: mm)

Model No.	Tubing ID	Material	ΦA	ΦB	ΦC	H
RH-G016	1.0~1.6	PP	1.60	0.99	2.80	18.6
RH-G024	1.6~2.4	PP	2.40	1.58	3.40	18.6
RH-G032	2.4~3.2	PP	3.20	2.46	4.20	19.6
RH-G040	3.2~4.0	PP	4.00	2.96	6.00	25
RH-G048	4.0~4.8	PP	4.80	3.94	6.80	25
RH-G064	4.8~6.4	PP	6.35	5.42	8.40	27.1

Model No.	Tubing ID	Material	ΦA	ΦB	ΦC	H
RH-M016	1.0~1.6	PP	1.60	0.99	2.80	19
RH-M024	1.6~2.4	PP	2.40	1.58	3.40	20
RH-M032	2.4~3.2	PP	3.20	2.46	4.20	20
RH-M040	3.2~4.0	PP	4.00	2.46	6.00	27.5
RH-M048	4.0~4.8	PP	4.80	2.96	6.80	27.5
RH-M064	4.8~6.4	PP	6.35	4.43	8.40	27.5

Model No.	Tubing ID	Material	ΦA	ΦB	ΦC	H	Thread	Stator	Nut
RH-C-M016	1.0~1.6	PP	1.60	0.99	2.40	25	1/4-28UNF	C1/4-28	LM-1/4-28P
RH-C-M024	1.6~2.4	PP	2.40	1.58	3.40	25	1/4-28UNF	C1/4-28	LM-1/4-28P
RH-C-M032	2.4~3.2	PP	3.20	2.46	4.20	25	1/4-28UNF	C1/4-28	LM-1/4-28P
RH-C-M040	3.2~4.0	PP	4.00	2.96	6.00	32.5	M10	C10	LM-M10P
RH-C-M048	4.0~4.8	PP	4.80	3.94	6.80	32.5	M10	C10	LM-M10P
RH-C-M064	4.8~6.4	PP	6.35	5.42	8.40	32.5	M10	C10	LM-M10P



Dimension (unit: mm)

