


MasterFlex Viton Tube



Viton is strong, resilient and flexible, provides excellent resistance to corrosive liquid such as acids, alkalis and some aggressive chemical solvents

- Continuous operating temperature: -32°C to +204°C
- Hardness: 75 A
- Ultra pure, non-toxic, excellent chemical inertness
- Breathability: low breathability
- UV resistance: A
- Ozone resistance: A
- Resistance to acids and bases: strong acids and bases (e.g. 98% concentrated sulphuric acid), some corrosive organic solvents (e.g. chloroform)
- Pump tube service life: 500 hours
- Opaque black color helps protect light sensitive fluid
- Certification: FDA and ADCF compliant

Rubber Tube / Viton Tube
Tubing

53
+
54

Model No.	ID	OD	WT (mm)	Colour
MasterFlex Viton Tube				
9763228 (1.02*2.62)	1.02	2.62	0.8	Opaque black
9763238 (1.65*3.25)	1.65	3.25	0.8	Opaque black
9763242 (2.06*3.66)	2.06	3.66	0.8	Opaque black
9763246 (2.54*4.14)	2.54	4.14	0.8	Opaque black
96412-13(13#)	0.8	4.00	1.6	Opaque black
96412-14(14#)	1.6	4.80	1.6	Opaque black
96412-16(16#)	3.2	6.40	1.6	Opaque black
96412-25(25#)	4.8	8.00	1.6	Opaque black
96412-17(17#)	6.4	9.60	1.6	Opaque black
96412-18(18#)	7.9	11.10	1.6	Opaque black