


PEEK One-Piece Fittings

Features:


- PEEK material, which insoluble in almost all solvents except concentrated sulfuric acid
- Superior to PTFE in terms of high temperature and abrasion resistance, high coefficient of friction and impact resistance
- High mechanical strength, fittings are not easily loosened after tightening
- Continuous operating temperature: -100°C~+250°C
- Original one-piece hand tightening fittings
- 10-32 threaded coned head for HPLC system application
- Suitable for 1.6 mm OD PEEK/PTFE/ETFE/FEP and other semi-rigid tubes

★ Application Note




Nut


Unit: mm




WK11A-10-32-K



WK11E-10-32-K



WK11C-10-32-K



WK11G-10-32-K

Model No.	Thread	Tube OD	Pressure Rating	Material
WK11A-10-32-K	10-32UNF Coned	1.6mm OD	8.64Mpa	PEEK natural
WK11E-10-32-K	10-32UNF Coned	1.6mm OD	8.64Mpa	PEEK natural
WK11C-10-32-K	10-32UNF Coned	1.6mm OD	8.64Mpa	PEEK natural
WK11G-10-32-K	10-32UNF Coned	1.6mm OD	8.64Mpa	PEEK natural