

## Saint-Gobain Rubber Tube



Saint-Gobain tubing is standard choice for peristaltic pumps

- Excellent mechanical strength, autoclavable, stretchable, non-toxic, nonhemolytic
- Continuous operating temperature: -60°C to + 135°C
- Hardness: 64 A
- Compatible with most CIP solutions and SIP
- Breathability: 80 (for nitrogen test, the higher the value, the better the tube breathability, and vice versa).
- UV resistance: A
- Ozone resistance: B
- Resistance to acids and alkalis: average
- Pump tube service life: 1000 hours (this value is provided by Saint-Gobain and is for reference only)
- Certification: USP Class VI, FDA, NSF, ISO 10993

| Model No.                  | ID   | OD    | WT (mm) | Colour       |
|----------------------------|------|-------|---------|--------------|
| Saint-Gobain Rubber Tube   |      |       |         |              |
| Pharmed BPT-0038-008       | 0.38 | 1.98  | 0.8     | Opaque beige |
| Pharmed BPT-0064-008       | 0.64 | 2.24  | 0.8     | Opaque beige |
| Pharmed BPT-0076-008       | 0.76 | 2.36  | 0.8     | Opaque beige |
| Pharmed BPT-0089-008       | 0.89 | 2.49  | 0.8     | Opaque beige |
| Pharmed BPT-0102-008       | 1.02 | 2.62  | 0.8     | Opaque beige |
| Pharmed BPT-0114-008       | 1.14 | 2.74  | 0.8     | Opaque beige |
| Pharmed BPT-0130-008       | 1.30 | 2.90  | 0.8     | Opaque beige |
| Pharmed BPT-0142-008       | 1.42 | 3.02  | 0.8     | Opaque beige |
| Pharmed BPT-0152-008       | 1.52 | 3.12  | 0.8     | Opaque beige |
| Pharmed BPT-0165-008       | 1.65 | 3.25  | 0.8     | Opaque beige |
| Pharmed BPT-0185-008       | 1.85 | 3.45  | 0.8     | Opaque beige |
| Pharmed BPT-0206-008       | 2.06 | 3.66  | 0.8     | Opaque beige |
| Pharmed BPT-0229-008       | 2.29 | 3.89  | 0.8     | Opaque beige |
| Pharmed BPT-0254-008       | 2.54 | 4.14  | 0.8     | Opaque beige |
| Pharmed BPT-0279-008       | 2.79 | 4.39  | 0.8     | Opaque beige |
| Pharmed BPT-0200-010       | 2.00 | 4.00  | 1.0     | Opaque beige |
| Pharmed BPT-0300-010       | 3.00 | 5.00  | 1.0     | Opaque beige |
| Pharmed BPT-0400-010       | 4.00 | 6.00  | 1.0     | Opaque beige |
| Pharmed BPT-0080-016 (13#) | 0.80 | 4.00  | 1.6     | Opaque beige |
| Pharmed BPT-0160-016 (14#) | 1.60 | 4.80  | 1.6     | Opaque beige |
| Pharmed BPT-0240-016 (19#) | 2.40 | 5.60  | 1.6     | Opaque beige |
| Pharmed BPT-0320-016 (16#) | 3.20 | 6.40  | 1.6     | Opaque beige |
| Pharmed BPT-0400-016       | 4.00 | 7.20  | 1.6     | Opaque beige |
| Pharmed BPT-0480-016 (25#) | 4.80 | 8.00  | 1.6     | Opaque beige |
| Pharmed BPT-0640-016 (17#) | 6.40 | 9.60  | 1.6     | Opaque beige |
| Pharmed BPT-0800-016 (18#) | 8.00 | 11.20 | 1.6     | Opaque beige |
| Pharmed BPT-0640-024 (24#) | 6.40 | 11.20 | 2.4     | Opaque beige |
| Pharmed BPT-0480-025 (15#) | 4.80 | 9.80  | 2.5     | Opaque beige |
| Pharmed BPT-0127-032 (82#) | 12.7 | 19.10 | 3.2     | Opaque beige |
| Pharmed BPT-0950-032 (73#) | 9.50 | 15.90 | 3.2     | Opaque beige |