

Silicone Tube

A translucent medical/food grade tubing which is odourless, non-toxic, Autoclavable, stretchable

- Continuous operating temperature -50°C to +200°C
- Hardness: 55 A
- Meets FDA standard, high sterility
- Anti-thermal oxidation, ozone Aging, light aging and weather aging
- Platinum cured for reduced residual levels of siloxanes and other volatiles
- Reduced spallation than rubber tube



Model No.	ID	OD	WT (mm)	Colour
Silicone Tube				
96405 (0.5*2.1)	0.50	2.10	0.8	Transparent
96406 (2.4*4)	2.40	4.00	0.8	Transparent
96410 (1*3)	1.00	3.00	1.0	Transparent
96411 (2*4)	2.00	4.00	1.0	Transparent
96413 (3*5)	3.00	5.00	1.0	Transparent
96421 (1.5*3)	1.50	3.00	1.5	Transparent
96402-13 (13#)	0.80	4.00	1.6	Transparent
96402-14 (14#)	1.60	4.80	1.6	Transparent
96402-19 (19#)	2.40	5.60	1.6	Transparent
96402-16 (16#)	3.20	6.40	1.6	Transparent
96402-25 (25#)	4.80	8.00	1.6	Transparent
96402-17 (17#)	6.40	9.60	1.6	Transparent
96402-18 (18#)	7.90	11.10	1.6	Transparent
96521 (1.6*6.4)	1.60	6.40	2.4	Transparent
96402-15 (15#)	4.80	9.80	2.5	Transparent
96402-24 (24#)	6.40	11.40	2.5	Transparent
96402-35 (35#)	7.90	12.90	2.5	Transparent
96402-36 (36#)	9.60	14.60	2.5	Transparent
16*22.4	16.00	22.40	3.2	Transparent
6.4*13	6.40	13.00	3.3	Transparent
96402-73 (73#)	9.60	16.60	3.5	Transparent
96402-82 (82#)	12.70	19.70	3.5	Transparent
BS-GJG-Φ10*Φ18	10.00	18.00	4.0	Transparent
12.7*22.3 (88#)	12.70	22.30	4.8	Transparent
25.4*35 (92#)	25.40	35.00	4.8	Transparent
9.5*22.1 (86#)	9.50	22.10	6.3	Transparent
19*31.6 (90#)	19.00	31.60	6.3	Transparent